




Solutions for Complex

Shoulder Girdle and Rib Injuries



**Clavicle Plating System
with Acu-Sinch® Repair System**
33 precontoured plating options,
plus a coracoclavicular anchor designed to
preserve the inferior cortex





Polarus® 3 Solution
The nail has a 4-degree curvature
designed to help preserve the rotator
cuff tendon and facilitate insertion

The plate has two bendable tabs
to support posterior comminution, and
calcar screws designed for medial
buttress support



Scapula Plating System
The industry's first fracture-specific
system, with comprehensive plating options
designed to streamline surgery

Glenoid Plate



Acute RibLoc® U+ Chest Wall
Plating System



For more information:
go.acumed.net/Shoulder
888.627.9957