


## From Straightforward to Specialized

## Our Hand Solutions Have You Covered

### OsteoMed Hand Plating System

A **full range of low-profile plates** in 4 size modules, hand fusion plates in 2 sizes, and innovative cannulated screws.



OsteoMed products



For more information:  
go.acumed.net/HS  
888.627.9957

JAD00-89-A | Effective: 2021/08 | © 2021 Acumed® LLC

### Acumed Hand Fracture System

**Specialized solutions** for challenging fractures and nearly 100 cut-to-fit plate options.

#### Customizable Standard Plates

Most plates can be cut to length and bent to fit a wide array of fractures. The included plate cutter is designed to create a smooth, rounded edge

#### 1.3 mm Rolando Fracture Hook Plate

A solution for a comminuted fracture of the base of the first metacarpal

#### Acutrak 2® Screws

The first and only continuously variable thread pitch headless compression screws on the market

#### 1.3 mm Metacarpal Neck Plate

The plate's 3 distally pointing converging screws support metacarpal head fixation

#### Small Bone External Fixator

Aids in reduction and compression to help realign fragmented bone

Acumed products